

Patient Study ID:

Visit number (circle): 0 2 5 8 12 24 Ad-hoc 1 Ad-hoc 2

ReceivedDate: Time: By: **Serum and urine separation****Tempus Blood RNA tube**

	Serum	Urine	Blood draw time:
Centrifugation start time:	<input type="text"/>	<input type="text"/>	<input type="text"/>
No. of aliquots frozen (max. 6):	<input type="text"/>	<input type="text"/>	Freezing time: <input type="text"/> (3 hours from blood draw)
Vol. of aliquots (ml):	0.5	0.9	
Freezing time:	<input type="text"/>	<input type="text"/>	

PBMC preparationVol. EDTA blood (ml): **Cell counts** (refer to local counting procedure)Cells/ml° = no. cells counted in 1mm² x dilution factor x 10⁴

° assumes depth of 0.1mm

Total no. of cells = cells/ml x total volume

Counting method used:Burker ☐Neubauer ☐Other (state)

	PBMC*	Monocytes	CD4 ⁺ T cells
Total cell vol.:	<input type="text"/>	<input type="text"/>	<input type="text"/>
No. cells counted:	<input type="text"/>	<input type="text"/>	<input type="text"/>
Dilution factor:	<input type="text"/>	<input type="text"/>	<input type="text"/>
Cells/ml:	<input type="text"/>	<input type="text"/>	<input type="text"/>
Total cell no.:	<input type="text"/>	<input type="text"/>	<input type="text"/>

Freezing cells

	PBMC (40% total)*	Monocytes	CD4 ⁺ T cells
Vol. to give 1x10 ⁵ cells (1x10 ⁵ ÷ cells/ml):	<input type="text"/> (white cap)	<input type="text"/> (amber cap)	<input type="text"/> (blue cap)
No. of aliquots frozen (orange caps):	6	1	1
No. of cells per aliquot (orange caps):	<input type="text"/>	<input type="text"/>	<input type="text"/>
Freezing time:	<input type="text"/>	<input type="text"/>	<input type="text"/>

Cells for flow cytometry purity check

	PBMC	Monocytes	CD4 ⁺ T cells
Vol. to give 2x10 ⁵ cells (2x10 ⁵ ÷ cells/ml):	<input type="text"/> x2	<input type="text"/> x1	<input type="text"/> x1

Monocyte and CD4⁺ T cell isolation**CD14 Microbeads**Lot no.: Expiry: **CD4 Microbeads**Lot no.: Expiry: **LS columns (x2)**Lot no.: Expiry:



Monocyte isolation

Total no. PBMC for cell isolation (60% total*):

No. PBMC for monocyte isolation:

(Subtract 4×10^5 cells used for the purity check
from total no. PBMC for cell isolation)

Vol. MACS buffer (80 μ l / 10×10^6 PBMC):
 μ l
Vol. CD14 MicroBeads (20 μ l / 10×10^6 PBMC):
 μ l
CD4⁺ T cell isolationAdded 160 μ l MACS Buffer? (Tick)
☐
Added 40 μ l CD4 MicroBeads? (Tick)
☐
CD4⁺ T cell lysate

QIAzol Lysis Reagent, Lot no.:

No. of lysed CD4⁺ T cells:

Freezing time:

Purity check by flow cytometry

Well/Tube	Cell type	Stain	Added 6 μ l antibody cocktail per well? (Tick)	Purity (%)
1	PBMC	Unstained	-	-
2	PBMC	Stained		-
3	Monocytes	Stained		
4	CD4 ⁺ T cells	Stained		

Antibodies to add to stained well:

Antibody	Vol. per test (μ l)	Vol. for 4x cocktail (μ l)	Lot no.
CD3 BV421	1 μ l	4 μ l	
CD56 FITC	1 μ l	4 μ l	
CD8 PE	1 μ l	4 μ l	
CD4 PE-Cy7	1 μ l	4 μ l	
CD19 APC	0.5 μ l	2 μ l	
CD14 APC-Cy7	1 μ l	4 μ l	

Notes:

Processed by:

Signature:

Date: